

SP200-1 Cell Magnetic Stirrer

Magnetic Series

Introduction

Cell Magnetic Stirrer is designed for research universities, laboratories and factories, used to culture and stir cells.

Features

1. Metal plate shell, firm and durable, thermostability and anti-corrosion.
2. TFT-LCD displays speed, starting time & end time, cycle numbers and programing pegments.
3. DC brushless motor makes stable operation, low noise, long using life.
4. Various working modes, including fixed speed stirring and cycle stirring.
5. Programmable. It can store and memorize nine frequently used settings S1~S9. One key to switch.
6. The control panel designed 30°slopes for convenient operation.
7. Magnetic stirring technology, stable running. Specially designed for cell culture with low speed.
8. Accelerate and decelerate gradually with low shear, keep cells growing safely.

Parameters

Model	SP200-1
Stirring Point Quantity	1
Max. Stir Capacity	3L
Speed Range	5~200rpm,step:1 rpm
Platform Material	stainless stee
Dimension	W.180 x D.283 x H.88mm
Voltage	AC 100~240V, 50/60Hz
Net weight	2.2kgs

